







## SUPERIOR EASE OF USE AND COMFORT

G+N anti-embolism stockings help to prevent the formation of Deep Vein Thrombosis (DVT) in the legs by applying clinically proven graduated compression.<sup>1,2</sup> G+N anti-embolism stockings are designed to be simple to fit and apply, which reduces sizing errors and nursing time.



DUOFIT<sup>™</sup> technology provides a two-way stretch for easier application.



Single ankle measurement reduces nursing time when measuring patients.



Six sizes within each length reduces inventory.



Colour coded packaging reduces time identifying correct size.







The compression within the garment

Anti-bacterial fabric for odour management.

Improved inspection hole

design and placement for


Available in below knee

superior comfort.

and thigh length.



### HOW DOES IT WORK?



### INTELLIGENT SOLUTIONS FOR DVT PREVENTION

With only 6 sizes within each style, G+N anti-embolism stockings allow hospitals to benefit from a reduced inventory - saving space in stockrooms and nursing time when identifying the correct size.

1. Independent Sigel Certification. Data on File.

2. Sigel B., et al. Type of Compression for Reducing Venous Stasis. Archives of Surgery. February 1975. Vol 110: 171-175.



# Only 6 sizes with colour coded packaging

### PRODUCT CODES & SIZING GUIDE



<b>FITLEGS</b> <sup>®</sup>	Suitability Check					
Ankle measurement	Size		V KNEE Order Code		IGH Order Code	Maximum Thigh Circumference
16 - 19cm	Slimline	EGD9095	403002			
19 - 23cm	Small	EGD8322	401000	EGD8326	401004	65cm
23 - 26cm	Medium	EGD8323	401001	EGD8327	401005	70cm
26 - 29cm	Large	EGD8324	401002	EGD8328	401006	78cm
29 - 33cm	Ex-Large	EGD8325	401003	EGD8329	401007	82cm
33 - 35cm	Bariatric	EGD8338	403000			

#### FITLEGS® AES GRIP

Suitability Check

Ankle measurement	Size		V KNEE Drder Code		IGH Drder Code	Maximum Thigh Circumference
16 - 19cm	Slimline	EGD9096	403003			
19 - 23cm	Small	EGD8330	402000	EGD8334	402004	65cm
23 - 26cm	Medium	EGD8331	402001	EGD8335	402005	70cm
26 - 29cm	Large	EGD8332	402002	EGD8336	402006	78cm
29 - 33cm	Ex-Large	EGD8333	402003	EGD8337	402007	82cm
33 - 35cm	Bariatric	EGD8339	403001			

Tape measure – Order code: 405001																															

(E

Griffiths and Nielsen Ltd Maydwell Avenue, off Stane Street, Slinfold, Horsham, West Sussex RH13 OGN Tel: +44 (0) 1403 799 190 Fax: +44 (0) 845 263 8907 medsales@gandn.com www.gandn.com

