



The original waterproof protector



Helping patients shower and bathe with confidence.

Long term protection for:

- Plaster casts
- Dressings
- Bandages
- PICC lines
- Compression hosiery

Versatile:

- Shower and bath
- One year guarantee
- Slip resistance
- Comprehensive sizing options

Available on:

- FP10 Prescription
- NHS Supply Chain
- SBS Framework



Designed and manufactured in the UK since 1994.

www.limboproducts.co.uk



"Since using the Limbo sleeves for our PICC patients we have had fewer PICCs falling out and our infection numbers have fallen."

Meinir Hughes - Intravenous Access Nurse Specialist

With Limbo your patients are totally covered:

Improve compliance and reduce costs to the NHS

- Reliable waterproof protection with a Limbo improves patient wellbeing, compliance and outcomes.
- Avoid wasted resources resulting from wet bandages, plaster casts, PICC lines or compression hosiery.

Easy to use and effective - totally covered

- Soft, flexible neoprene seal is easy to put on and sits flat against the skin providing a complete waterproof barrier.
- Non-constrictive and comfortable.
- Durable and versatile, a long-term solution for a wide variety of patient needs.

Shower AND bath - totally covered

- Not just showerproof - Limbos can be submerged, ideal for the bath.
- Perfect for hydrotherapy sessions.

Limbo's for all limb sizes - totally covered

- Simple size selector.
- Seal circumferences up to 90cm to help more of your patients.
- Customised service available for amputees or external fixator frames.

"Limbo has consistently proven to be a reliable and versatile solution. Their range of sizes means we can confidently support individuals of all shapes and conditions; its durability ensures repeated use without compromising effectiveness. Most importantly, Limbo helps maintain patients' dignity and independence, significantly improving their quality of life during recovery. An invaluable tool in both acute and community settings."

Luxmi Dhoonmoon
Nurse Consultant Tissue Viability
Queen's Nurse



One year guarantee - quality and value totally covered

- UK made with premium materials, individually inspected for quality.
- Reusable with an unrivalled one year guarantee - perfect for long term conditions.
- Increase compliance, avoid repeat prescriptions and associated costs.

Built in slip resistance - patient safety totally covered

- Triple layer sleeve purposely designed for better traction.
- Latex free, a trusted and safe solution since 1994.

Patient and clinical support team - totally covered

- West Sussex team available on phone and email for information and advice.
- Product demonstrations on educational days in person or virtually.
- Free samples available on request.
- Patients can also buy from our website or phone team. Rapid delivery with Royal Mail Tracked 24.

Prescribe and order via FP10 prescription, NHS Supply Chain and SBS Framework. Order codes and contact details on reverse.

"Limbo has a wonderful track record, exceptional patient support and fast delivery. Feedback is always positive, it's easy to use, reliable and cost effective - no more wet dressings, ruined casts and compromised wound sites, reducing the risk of infection - nothing else comes close! Limbo's educational support is exemplary too."

Maria A. Hughes
Nurse Consultant Tissue Viability
Queen's Nurse



NHS Supply Chain order codes and PIP codes

Product Description	Patient Weight	Measurements. Height, limb or cast circumference	MPC	NHS SC CODE	PIP CODE
Foot					
Adult Foot (slim-medium)	Max shoe size 11	Ankle 20-25cm (8-10"). Max stretch 33cm	M20	EAZ 723	416-6476
Adult Foot (medium-large)	Max shoe size 11	Ankle 23-34cm (9-13"). Max stretch 40cm	M25	EAZ 701	416-6484
Half Leg					
Adult Half Leg (slim)	Up to 64kg (10st)	165cm and above (5'5")	35cm - 39cm	M76	EAZ 710
Adult Half Leg (standard)	64-102kg (10-16st)		39cm - 54cm	M80	EAZ 718
Adult Half Leg (large)	102-140kg (16-22st)	Under 165cm (5'5")	52cm - 57cm	M180	EAZ 709
Adult Half Leg (slim-short)	Up to 64kg (10st)		35cm - 39cm	M76S	EAZ 702
Adult Half Leg (short)	64-102kg (10-16st)		39cm - 54cm	M80S	EAZ 707
Adult Half Leg (large-short)	102-140kg (16-22st)		52cm - 57cm	M180S	EAZ 717
Adult Half Leg (extra heavy)	140kg+ (22st+)	153-183cm (5'0"- 6'0")	65cm - 90cm	M200HL	EAZ 719
Full Leg					
Adult Full Leg (standard)	64-102kg (10-16st)	165-183cm (5'5"- 6'0")	52cm - 65cm	M100	EAZ 720
Adult Full Leg (standard short)	64 to 102kg (10-16st)	Under 165cm (5'5")	52cm - 65cm	M100S	N/A
Adult Full Leg (standard long)	64 to 102kg (10-16st)	Over 183cm (6'0")	52cm - 65cm	M100L	N/A
Adult Full Leg (slim)	Up to 64kg (10st)	under 183cm (under 6'0")	39cm - 54cm	M80LL	EAZ 722
Adult Full Leg (extra large)	102kg+ (16st+)	153-183cm (5'0"- 6'0")	65cm - 90cm	M200FL	EAZ 721

Product Description	Patient Weight	Measurements. Limb or cast circumference	MPC	NHS SC CODE	PIP CODE
Elbow / PICC					
Adult Elbow (extra slim)	Up to 51kg (8 stone)	22-25cm (8.5-10") bicep	M45	EAZ 85077	430-4960
Adult Elbow (slim-medium)	51-70kg (8-11 stone)	25-29cm (10-11.5") bicep	M65	EAZ 705	386-4303
Adult Elbow (medium - large)	70kg+ (11 stone +)	30-39cm (12-15") bicep	M75	EAZ 711	386-4311
Adult Elbow (Extra large)	108kg+ (17 stone +)	39-54cm (12-21") bicep	M85	EAZ 85076	430-4978
Hand					
Adult Hand	Slim to medium build	22-25cm (8.5-10") forearm	MITT	EAZ 85075	N/A
Adult Arm					
Adult Half Arm (slim-medium)	51-70kg (8 - 11 stone)	Limb 25-29cm, lower bicep, cast 40cm	M60	EAZ 706	N/A
Adult Half Arm (medium-large)	70kg + (11 stone+)	Limb 30-39cm, lower bicep, cast 44cm	M67	EAZ 699	N/A
Adult Full Arm (standard)	70-108kg (11-17 stone)	Limb 30-39cm, top of bicep, cast 44cm	M70	EAZ 715	431-9158
Adult Full Arm (slim)	51-70kg (8-11 stone)	Limb 25-29cm, top of bicep, cast 40cm	M60L	EAZ 716	431-9166
Adult Full Arm (large)	108kg+ (17 stone+)	Limb 39-54cm, top of bicep, cast 57cm	M87	EAZ 703	431-9141
Child					
Child Full Leg 6 to 7 Years	Size is typical for age; select as you would do for clothing.	Limb 35-39cm, cast 44cm, length 61cm	FL67	EAZ 714	N/A
Child Full Leg 8 to 10 Years		Limb 39-54cm, cast 57cm, length 73cm	FL810	EAZ 704	N/A
Child Full Leg 11 to 13 Years		Limb 39-54cm, cast 57cm, length 83cm	FL1113	EAZ 700	N/A
Child Full Arm 6 to 7 Years		Limb 22-25cm, cast 33cm, length 49cm	FA67	EAZ 712	N/A
Child Full Arm 8 to 10 Years		Limb 22-25cm, cast 33cm, length 58cm	FA810	EAZ 708	N/A
Child Full Arm 11 to 13 Years		Limb 25-34cm, cast 40cm, length 58cm	FA1113	EAZ 713	N/A

For more information or for patients wishing to purchase a Limbo, please contact us at:

www.limboproducts.co.uk Email: info@limboproducts.co.uk Phone: 01243 573417

Scan to visit our website for more information

LB01 SEPT25



The original waterproof protector

Unrivalled benefits and value for patients and the NHS.