



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant



Technical sheet



Skin Stapler

Reference(s)

F-35WB – wide staples
F-35RB – regular staples

Description

Disposable Skin Stapler with 35 stainless steel staples

Model	F-35RB	F-35WB
Staple diameter	0.5mm	0.6mm
Staple size after closure	5.9 mm*3.9mm 	7.2mm*4.9mm 
Qty of Staples (ea)	35	35
Sterilization method	EtO	EtO

Main materials

Body: Acrylonitrile butadiene styrene (ABS)
Staples cartridge: Polycarbonate
Staples: Medical-grade stainless steel

Pack factor

5 units per dispenser
100 units per carton (20 dispenser boxes)

Packaging

Peel packaging

Legal Manufacturer

Ningbo Advan Electrical Co.,Ltd
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: IIa sterile
• Please refer to complete labeling and instruction for use (IFU) before using the product.



STERILE EO

© 2022 Medline Industries LP or one of its subsidiaries. All rights reserved.
TDS_Skinstapler_EN04