



Peripheral neuropathy screening device

Pocket-sized **dependability**



OWEN MUMFORD

Neuropen[®]

Diabetic peripheral neuropathy is one of the most common complications of diabetes, potentially leading to loss of sensation, foot ulcerations and eventually amputation¹.

Treatment of these complications can be expensive² yet many are avoidable³, making early identification and intervention key to improving patient outcomes and reducing health costs.

Neuropen[®] is the only peripheral neuropathy screening device incorporating both pressure and sharpness tests, giving clinicians a simple, compact solution for testing large and small nerve fibres in the foot.

Neuropen[®]

The compact main body of the **Neuropen**[®] is designed to protect the monofilament between uses.



Neurotips[™]



40g pressure test

Neurotips[™] accurately assess reduced sensation to sharpness/pain in small nerve fibres. When used with **Neuropen**[®], **Neurotips**[™] ensure that a consistent force of 40g can be exerted safely onto the skin. **Neurotips**[™] are supplied irradiated and are for single use.



Accurate

Clinicians can be confident in achieving reliable results with **Neuropen® Monofilaments** which have been shown to be amongst the most accurate when compared to other manufacturers⁴.



Neuropen® Monofilaments

10g pressure test

Neuropen® Monofilaments accurately assess touch/pressure sensation in large nerve fibres, enabling clinicians to map areas of reduced pressure perception. The 10g monofilaments exert a specific, repeatable force on the test site.

They can be replaced without needing to discard the whole device.

Neuropen®

Dual function screening device with replaceable monofilaments



Neuropen®

1 x Neuropen; 1 x 10g Neuropen Monofilament; 1 x Neurotip

Starter pack:
OM Code: P005138
PIP Code: 270-6281



Neuropen® Monofilaments

10g Pressure test

Box of 5:
OM Code: P005136
PIP Code: 270-6299



Neurotips™

40g Pressure test

Box of 100:
OM Code: P005019
PIP Code: 270-6273

For more information please visit us at owenmumford.com

Our accreditations include:

- Medical Device Directive 93/42/EEC
- MDSAP (ISO 13485:2016)
- Environmental Management Systems ISO14001



Certified



1. Guttormsen K, Chadwick P (2017) Diabetic neuropathy: Beyond the basics. Journal of Diabetes Nursing 21: 17–22 2. Smith AG, Živković S. The hidden costs of painful diabetic neuropathy revealed. Neurol Clin Pract. 2020 Feb;10(1):3-4. doi: 10.1212/CPJ.0000000000000759. Epub 2019 Oct 29. PMID: 32185100; PMCID: PMC7057076 3. State of the nation 2016, Diabetes UK. 4. Differences in the performance of commercially available 10-g monofilaments. J Booth, M J Young. Diabetes Care Jul 2000, 23 (7) 984-988.

© Copyright Owen Mumford 2024
NEUROPEN-UK-LT-EN/MICC/0224/01