



PRODUCT INFORMATION

| | |
|----------------------------|------------------|
| Item Number | GN91 |
| Product Colour | Cool Blue |
| Size Range | XS-XL |
| Fingertip Thickness | 0.08 mm |
| Palm Thickness | 0.06 mm |
| Typical Weight | 3.20 g |
| Shelf Life | 5 Years |

FEATURES

Chemical Resistant: Achieves type C chemical resistance according to EN ISO 374-1, for chemical protection against at least 1 chemical

Virus/Non Virus Standard: Conforms with the requirements of EN ISO 374-5, for protection against bacteria, fungi and virus

Flexible material: For improved feel and fit reducing fatigue during extended wear

Latex Free: Reducing the risk of hypersensitivity and allergic reactions

Low in Residual Chemicals: To reduce the risk of contact dermatitis associated with these chemicals

Manufactured From: Specially formulated nitrile butadiene polymer

Powder Free: Reduces the risk of dust contamination

Textured Finish: For improved grip

SECTORS

Automotive • Construction & Utility • Engineering and Manufacturing • Janitorial • Medical, Vet and Dental • Oil, Gas & Renewables



STANDARDS

AQL 1.5


This product is classified as Class I Medical Device in compliance with the requirements of Medical Device Regulation (EU) 2017/745 & UK MDR 2002 (SI2002 No.618 as amended) (UKCA)

EU Certified

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European standards EN ISO 21420, EN ISO 374-1, EN ISO 374-4 & EN ISO 374-5



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555



CAUTION

Contains accelerator chemicals which may cause type IV skin irritation and contact dermatitis in sensitised individuals.

PACKING

200 gloves per dispenser
10 dispensers per case
50 cases per pallet

INNER PACKAGING



STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. You are advised to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear.

DISPOSAL INSTRUCTIONS

Dispose of responsibly and in line with local guidance.

OUTER PACKAGING



| Size | Inner Packaging Barcode |
|--------|-------------------------|
| XS | 5024951016045 |
| Small | 5024951016052 |
| Medium | 5024951016069 |
| Large | 5024951016076 |
| XL | 5024951016083 |

| Size | Outer Packaging Barcode |
|--------|-------------------------|
| XS | 5024951016090 |
| Small | 5024951016106 |
| Medium | 5024951016113 |
| Large | 5024951016120 |
| XL | 5024951016137 |

| Inner Packaging | |
|--------------------------|------------------|
| Inner dimensions (LxWxH) | 230 x 125 x 73mm |

| Outer Packaging | |
|--------------------------|-------------------|
| Outer Dimensions (LxWxH) | 375 x 262 x 240mm |

PLEASE NOTE:

Weight does not form part of the specification and is subject to change. Palm thicknesses is approximate and based on an agreed manufacturing tolerance. All thickness and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555