


Product Specifications Sheet

Specification Sheet No	SSN-FT1358-015		
Product Name	Surgical Drape		
Product Description	C-Section Drape Set (Teardrop Fenestrated)		
Range	365 Healthcare		
Manufacturer	Sino Protection (HeFei) Medical Products Co., Ltd.		
Material	SPI TB 3 layers material, PE with viscose reinforcement, SPI TB material, SM PAD, Clear PE, white Spun lace, Cotton		
Sterilisation	EO Sterile		
Issue Date	31 st July 2025	Version No	R02



Intended use	To be used to drape the patient to isolate a site from other areas of the patient's body or to contribute to reducing the risk of cross contamination.
---------------------	--

Packaging					
Product code	Item	Dimension (cm)	Colour	Items Per Pack	Qty – Shipper Carton
84-006365	1 x C-Section Drapes	257/196 x 325 (±2)	Blue & Clear	1	5
	2 x Tape Strips	2 x 30			
	4 x Hand Towels	35.5 x 67.5			
	1 x Table Cover	140 x 205			
	1 x Baby Blanket	90 x 110	White		
	1 x Baby Hat	12 x 15			

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	Medical Devices Directive (EU) 94/32/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	No known		
Markings	365 Healthcare, Ref, Quantity, Product name, Product list, Image & dimension, Barcode, Manufacturing & Expiry date, Lot, CE Mark, Single use, Latex free, Medical device, EO Sterile, keep dry, keep away from sunlight, Do not use if package is damaged, EC-REP, Manufacturer name & address		
Additional information	This product is latex free.		
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	24/07/2024	1 st version	Jack Wong
R02	31/07/2025	Update baby hat dimension	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.