


Product Specifications Sheet

| | | | | |
|-------------------------------|--------------------------------------|-------------------|-----|---|
| Specification Sheet No | SSN-JFRU06-007 | | |  |
| Product Name | Suction | | | |
| Product Description | UHS Yankauer Sucker - Midi No Vacuum | | | |
| Range | Universal | | | |
| Manufacturer | Romsons International | | | |
| Material | Poly Vinyl Chloride | | | |
| Sterilisation | EO Sterile | | | |
| Issue Date | 2 nd September 2025 | Version No | R03 | |

| | |
|---------------------|--|
| Intended use | Intended for preoperative secretions, blood and surgical debris removal. |
|---------------------|--|

| Packaging | | | | |
|---------------|----------------|--------------------------|----------------|----------------------|
| Product code | Length (±10mm) | Outer Dimension (±0.5mm) | Items Per Pack | Qty – Shipper Carton |
| UNY410-100 | 250cm | Proximal end 8mm | 1 | 100 |
| UNY410 VC-100 | 250cm | Distal end 4mm | 1 | 100 |
| UNY411-100 | 250cm | Bubble diameter 12.5mm | 1 | 100 |
| UNY411 VC-100 | 250cm | Proximal end 8mm | 1 | 100 |
| UNY412-100 | 250cm | Distal end 4.8mm | 1 | 100 |
| UNY412 VC-100 | 250cm | Bubble diameter 12.5mm | 1 | 100 |
| UNY413-100 | 250cm | Proximal end 8mm | 1 | 100 |
| UNY413 VC-100 | 250cm | Distal end 6mm | 1 | 100 |
| | | Bubble diameter 12.5mm | 1 | 100 |

| | | | |
|-----------------------------------|--|-------------|---|
| Quality Certification | ISO 13485:2016 | | |
| Device Classification | Class IIa | Rule | 6 |
| Conformity Route | Medical Devices Directive (93/42/EEC) | | |
| Storage conditions | Store in a cool dry place away from heat and direct sunlight | | |
| Shelf life | 5 Years from manufacturing date if stored as recommended | | |
| Contraindications/Warnings | To be used by a qualified / a well-trained medical person. | | |
| Markings | Ref, Product name, Lot, Quantity, Barcodes, Manufacture date, Expiry date, CE Mark, Single use, Do not re-sterilise, Latex free, EO-sterile, DEHP free, Medical device, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer Name & Address, Distributor contact, Non-pyrogenic | | |
| Additional information | Product is packed in a peelable blister pack per unit. Product is Latex-free | | |

Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

| | |
|--|----------------------------------|
| 365 Healthcare | Mediq House |
| T: +44 (0) 1332 413 700 | Arundel Avenue, Castle Donington |
| W: https://mediq.co.uk/ | Derby, DE74 2HJ |

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

| Document Reference | Document Name | Reviewed By | Approved By | Revision Number | Issue Date |
|--------------------|-----------------------------|-------------|--------------|-----------------|-----------------------------|
| DMS 04-004 | Product Specification Sheet | Jack Wong | Racheal Begg | R02 | 03 rd April 2024 |

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

| Revision Number | Revision Date | Reason for Revision | Updated by |
|-----------------|---------------|--|------------|
| R01 | 29/02/2024 | Update to Mediq Healthcare UK brand replacing DSN | Jack Wong |
| R02 | 21/05/2024 | Update Mediq Healthcare UK address | Jack Wong |
| R03 | 02/09/2025 | Update with other product code within the same range | Jack Wong |
| | | | |
| | | | |

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

| Document Reference | Document Name | Reviewed By | Approved By | Revision Number | Issue Date |
|--------------------|-----------------------------|-------------|--------------|-----------------|-----------------------------|
| DMS 04-004 | Product Specification Sheet | Jack Wong | Racheal Begg | R02 | 03 rd April 2024 |