

**Multi Cup  
116AA200  
NPC Code FSX 014  
Quantity 200 per Carton**




KM 650633



**Tested to Comply with  
PAS29:1999**


**Outer Case Dimensions 360 x 430 x 290mm**



(01)00800239011163

**OUTER BAR CODE**

- SINGLE USE
- MOULDED FIBRE PULP
- MACERATABLE
- MADE FROM RECYCLED PRODUCT

- ❓ Ideal for taking urine/stool/sputum samples; dispensing tablets or as a measuring scoop (not for use with detergents)
- ❓ Disposable single use product to reduce risk of cross infection.
- ❓ Moulded products made from 100% recycled cellulose fibres.
- ❓ Fully disposable in Vernacare Disposal Units.

**Regulatory Compliance**

- ❓ Compliant with: Medical Device Regulation 2017/745 Class 1 (non-sterile)
- ❓ Product Specification Assessment PAS 29:1999
- ❓ Disposability: By Maceration
- ❓ ISO 9001:2015 / ISO 14001:2015

**Description:**

- ❓ Capacity: 200ml (max)
- ❓ Material: 100% Paper Pulp
- ❓ Dimensions: 160 x 90 x 65mm
- ❓ Weight: 8.1g (+/- 1g)
- ❓ Colour: Grey

**Fluid Retention Testing**

- ✓ Tested in Accordance with PAS 29:1999
- ✓ Guaranteed to hold fluid for 4 Hours
- ✓ Resistance Temperature: 35 °C (+/- 3)



**Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.**

+44 1204 529494 (UK) | (416) 661 5552 (N.America) | (800) 268 2422 (Toll Free)  
Vernacare Ltd, Matrix Park, 1 Western Ave, Chorley Lancs UK, PR7 7NB